## TECHNICAL DATA SHEET



## ANTIQUAX PRODUCT NAME Silver Bath

Waterloo Road, London NW2 7TX Telephone: +44(0)20 8450 4666

## CODE ANTQSBBX

Туре	Sachet.
Intended Use	To remove tarnish from solid silver and silver plate items.
Special Features	No rubbing required. Cleans and shines in minutes. Suitable for engraved and filigree pieces.  Non toxic
Preparation	Place sachet in a large plastic container. Place items to be cleaned in the container making sure that the items are touching each other and one item is in contact with the silver side of the sachet. Add hot water so as to cover the items and leave for 5 minutes. Remove items and wash in soapy water and dry. The items will now be free of all tarnish.
Application	A single application only required.
Coats	N/A
Drying Time	N/A
Re-coating Time	N/A
Application Temp.	Hot water temperature.
Coverage	N/A
Maintenance	Re-treat when necessary
Compatibility	Do not use on lacquered or artificially oxidised silver ware, candlesticks with weighted bases, knives with weighted handles, candlesticks with weighted bases, or articles with porous attachments such as wood, leather, horn, enamel, bone or baize.
Settlement	N/A
Shelf-life/Storage	Store in a cool dry place.
Hazards	Avoid contact with eyes skin. Persons with sensitive skin should wear rubber gloves. Keep away from children.
Clean-Up	Warm soapy water.
Sizes	3 x 50 gm sachets per box.